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Telephone:03-5220-8953

Facsimile :03-5220-8956

Endospecey[®] ES-50M Kit
Toxicolor[®] LS-50M Kit
PyroColor Diazo Reagents DIA 150-MP

Endotoxin test is based on the activation and gelation of LAL (Limulus Amebocyte Lysate), an extract of circulating amebocytes of horseshoe crabs (*Limulus polyphemus*, *Tachypleus tridentatus*, etc.), by endotoxin (Et) derived from Gram-negative bacteria. Since this reaction occurs as a series of enzymatic reactions, a chromogenic assay has been developed which utilizes as the index a color that is developed as a result of the enzymatic hydrolysis of a synthetic substrate.

Endospecey[®] ES-50M Kit and Toxicolor[®] LS-50M Kit are reagent kits for chromogenic assay of endotoxin. Each kit consists of a LAL reagent comprising lysate of *Limulus polyphemus* and a synthetic chromogenic substrate (Boc-Leu-Gly-Arg-pNA), a buffer solution for dissolving LAL reagent and distilled water (Et and β -glucan free).

Endospecey[®] ES-50M Kit is formulated so that it reacts specifically only with endotoxin without reacting with false-positive substances [substances having (1 \rightarrow 3)- β -D-glucan structure, such as fungal polysaccharides and washings of cellulose-type filter membranes].

Toxicolor[®] LS-50M Kit contains both factor C system (which reacts with endotoxin) and factor G system (which reacts with (1 \rightarrow 3)- β -D-glucan). If it is known that the sample does not contain false-positive substances with (1 \rightarrow 3)- β -D-glucan structure*¹, or if presence or absence of pseudo-positive substances is not of concern, use of Toxicolor[®] LS-50M Kit is recommended.

*¹ For the measurement of (1 \rightarrow 3)- β -D-glucan, we recommend the use of our product, Glucate[®] kit.

Characteristic features

- Conforms to Endotoxin Test of Japanese Pharmacopoeia.
- Offers a highly sensitive endotoxin assay.
- Applicable both to kinetic chromogenic assay using microplate and endpoint chromogenic assay.
- Endospecey[®] ES-50M Kit reacts specifically only with endotoxin.
- Supplied in multi packages convenient for multiple sample assay.

Lab wares required

Use the following recommended lab wares or those confirmed to be equivalent in material and function.

Et and β -Glucan free lab wares

Lab wares	Recommended lab wares	Code No.
Syringe	Toxipet Dispenser Syringe	900520
Tips	Toxipet Tip 200	900540
	Toxipet Tip 1000	900545
Microplate	Toxipet Plate LP	900570
Aluminum foil	Dry-heat sterilized aluminum foil	800804

Other lab wares required

Lab wares	Recommended lab wares	Code No.
Microplate reader	Wellreader MP-96	900625
Hot plate* ²		
Plate mixer* ²		
Test tube mixer	Test tube Mixer TM-251	900675

*² Not required if Wellreader MP-96 is used.

Components of each kit

020150 Endospecey [®] ES-50M Kit* ³	for 50 assays
③ Buffer	2.8 mL \times 1 vial
④ LAL reagent (factor C-pathway)	1 vial
⑤ Distilled water (Et and β -glucan free)	2.8 mL \times 1 vial

*³ Endotoxin specific

020130 Toxicolor [®] LS-50M Kit	for 50 assays
③ Buffer	2.8 mL \times 1 vial
④ LAL reagent	1 vial
⑤ Distilled water (Et and β -glucan free)	2.8 mL \times 1 vial

010145 PyroColor Diazo Reagents DIA150-MP	for 150 assays
(for endpoint colorimetry in microplates)	
(1A) HCl solution(for dissolving Reagent (1))	4.0 mL \times 2 vials
(1) Sodium nitrite	1.65 mg \times 2 vials
(2) Ammonium sulfamate (ammonium amidosulfate)	12 mg \times 2 vials
(3) <i>N</i> -(1-Naphthyl)ethylenediamine dihydrochloride	2.8 mg \times 2 vials
(3A) <i>N</i> -Methyl-2pyrrolidone solution (for dissolving Reagent (3))	4.0 mL \times 2 vials

Storage, expiry date

Storage: to be stored protected from light at 2 to 8°C (do not freeze)

Expiry date: indicated on the outer package

Options

- 010020 Control Standard Endotoxin (for Pyrochrome)
 010290 LRW (LAL Reagent Water)

Reconstitution of reagents

Reconstitution of LAL reagent:

Add the total volume (2.8 mL) of the buffer solution ③ to the vial of LAL reagent ④ of Endospecey[®] ES-50M Kit or Toxicolor[®] LS-50M Kit using a syringe (Et and β -glucan free), etc., cover the top of the vial with a sheet of aluminum foil (dry-heat sterilized), and mix and dissolve the content while avoiding foaming (do not use a test tube mixer since it causes foaming). Prepare LAL reagent immediately before use. Make sure that the reagent has dissolved completely. (Usually, it takes about 5 minutes to dissolve the reagent. While turbidity disappears and a clear solution is obtained.) LAL reagent should be unsealed and reconstituted before each use, and the reconstituted reagent should be for single use only.

Preparation and storage of standard solution:

For endotoxin tests conforming to Japanese Pharmacopoeia, use standard endotoxin (JP grade) as the standard and follow the procedures stipulated in the Pharmacopoeia for preparation and storage of standard solution.

Reconstitution of control standard solution:

Add 5.0 mL of LRW to the vial of Control Standard Endotoxin (for Pyrochrome) using a Toxipet tip 1000 (Et and β -glucan free) etc., cover the top of the vial with a sheet of aluminum foil (dry-heat sterilized), and mix the content vigorously using a test tube mixer for at least 1 minute while taking care so that the solution does not come contact with the aluminum foil. Then, make dilutions appropriately with LRW. Mix the content vigorously again immediately before use.

Concentration of the control standard solution:

Certificate of Analysis sheet indicating potency (EU/ng) of each lot of Control Standard Endotoxin standardized against USP-Reference Standard Endotoxin with each lot of Endospey ES-50M or Toxicolor LS-50M can be obtained on our website (<http://www.lalbiz.com/>). Based on the potency, calculate the endotoxin concentration in the reconstituted solution, and then make dilutions appropriately.

Storage of control standard solution:

For storing the standard solution, cover the top of the vial with a sheet of aluminum foil (dry-heat sterilized), seal the circumference with Parafilm, and store at 2 to 8°C. The standard solution is stable for 1 week.

Blank:

Use distilled water ⑤ (Et and β -glucan free) as the blank.

Preparation of diazo-coupling reagent:

Using the reagents supplied with the kit:

1. Add the total volume (4.0 mL) of HCl solution (1A) to the vial of sodium nitrite (1).
2. Reconstitute the vial of ammonium sulfamate (2) with 4.0 mL distilled water.
3. Add the total volume (4.0 mL) of *N*-Methyl-2-pyrrolidone solution (3A) to the vial of *N*-(1-Naphthyl)ethylenediamine dihydrochloride (3).
4. Mix and dissolve the contents by gently swirling the vials to prepare solutions (1), (2), and (3).

Blank value and assay range in standard assay procedure (The standard assay is carried out at 37°C for 30 minutes)

Blank value

Kinetic chromogenic: Change rate ≤ 1.00 mAbs/min

Endpoint chromogenic: Absorbance ≤ 0.100 Abs

Detection limit

< 0.001 EU/mL

Assay range (measurement by arithmetic scale)


Kinetic chromogenic: 0.002 - 0.1 EU/mL

Endpoint chromogenic: 0.002 - 0.1 EU/mL

Standard assay procedure


I. Kinetic chromogenic assay using a microplate

1. Transfer 50 μ L each of samples, standard solutions (serial dilutions) and the blank to wells of a microplate (Et and β -glucan free) using pipette tips (Et and β -glucan free).
2. Add 50 μ L of the reconstituted LAL reagent to each well containing the sample, standard solution or blank using a syringe (Et and β -glucan free).
3. Cover the microplate and place in a dedicated microplate reader. After stirring, measurement starts at 37°C [reading wavelength 405 nm, reference wavelength 492 nm] and Et concentrations in samples are calculated automatically after the reaction is completed. The standard assay reaction takes 30 minutes to complete.

II. Endpoint chromogenic assay using a microplate

1. Transfer 50 μ L each of samples, standard solutions (serial dilutions of Control Standard Endotoxin solution) and the blank to wells of a microplate (Et and β -glucan free) using pipette tips (Et and β -glucan free).
2. Add 50 μ L of the reconstituted LAL reagent to each well containing the sample, standard solution or blank using a syringe (Et and β -glucan free).
3. Cover the microplate, mix on a plate mixer, and incubate on a hot plate at 37°C for 30 minutes.
4. After assay incubation is completed, immediately add 50 μ L of Solution (1) to each well and gently shake the plate to mix. Promptly add solution (1) to all wells to ensure that all samples incubate for the same period of time. Add 50 μ L of Solution (2) to each well and shake plate to mix. Add 50 μ L of Solution (3) to each well and shake plate to mix. Full color (magenta) should develop immediately.
5. Place the microplate in a dual wavelength microplate reader and measure absorbance at 545 nm (reference wavelength 630 nm).

Precautions in use and handling

Please read this package insert carefully before using the kit.

1. All reagents contained in this kit are specified only for use in research and testing. It is not permitted to use them as pharmaceutical products, diagnostics or for any other purposes.
2. Use the kit according to the procedure described in the package insert. Use of the kit with procedures or for purposes other than described here cannot be guaranteed.
3. Do not use the reagent after its validity date has expired.
4. Since diazo-coupling reagent contains HCl solution and sodium nitrite, care should be exercised in handling the kit. Sodium nitrite is designated as an oxidative substance by Fire Service Law and other regulations. Material Safety Data Sheet (MSDS) is available upon request.
5. If the reagent gets into eyes or mouth, get first aid such as washing off the reagent thoroughly and get medical help if necessary.
6. After opening the vial, use the reagent as soon as possible. When storing, seal the vial and store according to the method specified. The performance of stored reconstituted LAL reagent cannot be guaranteed.
7. Use forceps to remove the aluminum cap of vial.
8. Alternate use of containers and accessories is not guaranteed.

SEIKAGAKU CORPORATION
6-1, Marunouchi 1-chome, Chiyoda-ku
Tokyo, Japan