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Facsimile :03-5220-8956

Endospecy® ES-24S Kit

The endotoxin test is based on the activation and gelation of LAL (Limulus Amebocyte Lysate), an extract of circulating amebocytes of horseshoe crabs (*Limulus polyphemus*, *Tachypleus tridentatus*, etc.), by endotoxin (Et) derived from Gram-negative bacteria. Since this reaction occurs as a series of enzymatic reactions, a chromogenic assay has been developed which utilizes as the index a color that is developed as a result of the enzymatic hydrolysis of a synthetic substrate.

Endospecy® ES-24S Kit is reagent kit for chromogenic assay of endotoxin. The kit consists of a LAL reagent comprising lysate of *Limulus polyphemus* and a synthetic chromogenic substrate (Boc-Leu-Gly-Arg-pNA), a buffer solution for dissolving LAL reagent and distilled water (Et and β -glucan free).

Endospecy® ES-24S Kit is formulated so that it reacts specifically only with endotoxin without reacting with false-positive substances [substances having (1 \rightarrow 3)- β -D-glucan structure, such as fungal polysaccharides and washings of cellulose-type filter membranes].*¹

*¹ For the measurement of (1 \rightarrow 3)- β -D-glucan, we recommend the use of our product, GlucateLL® kit.

Characteristic features

- Conforms to Endotoxin Test of Japanese Pharmacopoeia.
- Offers a highly sensitive endotoxin assay.
- Applicable both to single-type kinetic chromogenic assay using the vials as the reaction test tube and to endpoint chromogenic assay.
- Endospecy® ES-24S Kit reacts specifically with endotoxin only.
- Supplied in single packages permitting the use of the minimum amount of reagent required thus eliminating wastes.

Lab wares required

Et and β -Glucan free lab wares

Lab wares	Recommended lab wares	Code No.
Tips	Toxipet Tip 200	900540
	Toxipet Tip 1000	900545
Aluminum foil	Dry-heat sterilized aluminum foil	800804
Aluminum caps	Dry-heat sterilized aluminum caps	800807

Other lab wares required

Lab wares	Recommended lab wares	Code No.
Colorimeter	EG Reader SV-12* ² (for kinetic chromogenic assay)	900630
Test tube mixer	Test tube Mixer TM-251	900675

*² If EG Reader SV-12 is not used (endpoint chromogenic assay), a general purpose colorimeter, a thermostat, bath and an ice-water bath are required.

Components of kit

020170 Endospecy® ES-24S Kit	for 24 assays
③ Buffer	2.8 mL \times 4 vials
④ LAL reagent (factor C-pathway)	24 vials
⑤ Distilled water (Et and β -glucan free)	2.8 mL \times 2 vials

Storage, expiry date

Storage: to be stored protected from light at 2 to 8°C (do not freeze)

Expiry date: indicated on the outer package

Options

- 010020 Control Standard Endotoxin (for Pyrochrome)
- 010290 LRW (LAL Reagent Water)
- 010045 PyroColor Diazo Reagents DIA60-STV (required for endpoint chromogenic assay but not for kinetic chromogenic assay)

Reconstitution of reagents

Reconstitution of LAL reagent:

Add 200 μ L of the buffer solution ③ to the vial of LAL reagent ④ contained in Endospecy® ES-24S Kit using a syringe (Et and β -glucan free), etc., cover the top of the vial with an aluminum cap (dry-heat sterilized), and mix about 2 seconds with a test tube mixer and dissolve the content while avoiding foaming. Prepare LAL reagent immediately before use. Make sure that the reagent has dissolved completely. LAL reagent should be unsealed and reconstituted before each use, and the reconstituted reagent should be for single use only.

Preparation and storage of standard solution:

For endotoxin tests conforming to Japanese Pharmacopoeia, use standard endotoxin (JP grade) as the standard and follow the procedures stipulated in the Pharmacopoeia for preparation and storage of standard solution.

Reconstitution of control standard solution:

Add 5.0 mL of LRW to the vial of Control Standard Endotoxin (for Pyrochrome) using a Toxipet tip 1000 (Et and β -glucan free) etc., cover the top of the vial with a sheet of aluminum foil (dry-heat sterilized), and mix the content vigorously using a test tube mixer for at least 1 minute while taking care so that the solution does not come contact with the aluminum foil. Then, make dilutions appropriately with LRW. Mix the content vigorously again immediately before use.

Concentration of the control standard solution:

Certificate of Analysis sheet indicating potency (EU/ng) of each lot of Control Standard Endotoxin standardized against USP-Reference Standard Endotoxin with each lot of Endospecy ES-50M or Toxicolor LS-50M can be obtained on our website (<http://www.lalbiz.com/>). Based on the potency, calculate the endotoxin concentration in the reconstituted solution, and then make dilutions appropriately.

Storage of control standard solution:

For storing the standard solution, cover the top of the vial with a sheet of aluminum foil (dry-heat sterilized), seal the circumference with Parafilm, and store at 2 to 8°C. The standard solution is stable for 1 week.

Blank:

Use distilled water ⑤ (Et and β -D-glucan free) as the blank.

Preparation of diazo-coupling reagent:

Using the reagents supplied with the kit:

1. Add the total volume (12.0 mL) of HCl solution (1A) to the vial of sodium nitrite (1).
2. Reconstitute the vial of ammonium sulfamate (2) with 12.0 mL distilled water.
3. Reconstitute the vial of *N*-(1-Naphthyl)ethylenediamine dihydrochloride (3) with 12.0 mL of distilled water.
4. Mix and dissolve the contents by gently swirling the vials to prepare solutions (1), (2), and (3).

Blank value and assay range in standard assay procedure

(The standard assay is carried out at 37°C for 30 minutes)

Blank value

Kinetic chromogenic: Change rate ≤ 3.00 mAbs/min

Endpoint chromogenic: Absorbance ≤ 0.200 Abs

Detection limit

< 0.001 EU/mL

Assay range (measurement by arithmetic scale)


Kinetic chromogenic: 0.002 - 0.1 EU/mL

Endpoint chromogenic: 0.002 - 0.1 EU/mL

Standard assay procedure

I. Kinetic chromogenic assay using EG Reader SV-12

1. Prepare the required number of vials containing reconstituted LAL reagent.
2. Transfer 200 μ L each of samples, standard solutions (serial dilutions) and the blank to vials containing LAL reagent using pipette tips (Et and β -glucan free) and stir on a test tube mixer about 2 seconds.
3. Place the LAL reagent vials in an EG Reader SV-12. Measurement starts at 37°C [reading wavelength 405 nm, reference wavelength 492 nm] and Et concentrations in samples are calculated automatically after the reaction is completed. The standard assay reaction takes 30 minutes to complete.

Precautions in use and handling

Please read this package insert carefully before using the kit.

1. All reagents contained in this kit are specified only for use in research and testing. It is not permitted to use them as pharmaceutical products, diagnostics or for any other purposes.
2. Use the kit according to the procedure described in the package insert. Use of the kit with procedures or for purposes other than described here cannot be guaranteed.
3. Do not use the reagent after its validity date has expired.
4. Since diazo-coupling reagent contains HCl solution and sodium nitrite, care should be exercised in handling the kit. Sodium nitrite is designated as an oxidative substance by Fire Service Law and other regulations. Material Safety Data Sheet (MSDS) is available upon request.
5. If the reagent gets into eyes or mouth, get first aid such as washing off the reagent thoroughly and get medical help if necessary.
6. After opening the vial, use the reagent as soon as possible. When storing, seal the vial and store according to the method specified. The performance of stored reconstituted LAL reagent cannot be guaranteed.
7. Use forceps to remove the aluminum cap of vial.
8. Alternate use of containers and accessories is not guaranteed.

SEIKAGAKU CORPORATION
6-1, Marunouchi 1-chome, Chiyoda-ku
Tokyo, Japan

II. Endpoint chromogenic assay using general purpose colorimeter

1. Prepare the required number of vials containing reconstituted LAL reagent and immerse them in an ice-water bath.
2. Transfer 200 μ L each of samples, standard solutions (serial dilutions) and the blank into vials containing reconstituted LAL reagent using pipette tips (Et and β -glucan free) and stir on a test tube mixer about 2 seconds.
3. Incubate LAL reagent vials in a thermostat bath at 37°C for 30 minutes.
4. After assay incubation is completed, immediately add 500 μ L of Solution (1) to each tube and vortex 2 seconds to mix. Promptly add solution (1) to all tubes to ensure that the samples reactions incubate for the same period of time. Add 500 μ L of Solution (2) to each tube and vortex 2 seconds to mix. Add 500 μ L of Solution (3) to each tube and vortex 2 seconds to mix. Full color (magenta) should develop immediately. The procedures need not be Et and β -glucan free. Measure absorbance at 545 nm using distilled water as the reference.